

Ortho Supply List

(Available through Henry Scheine, Patterson)

For Bonding Attachments:

Liquid Etch 23g Roughens enamel to create a bondable surface for attachments

Assure Plus Adhesion primer that improves composite retention on enamel

Transbond Supreme LV Flowable | 2 pk Composite used to build and bond attachments

For IPR:

Mosquito Diamond Tip IPR Burs Removes small amounts of enamel between teeth to create space for movement

IPR Strips 8mm Manual alternative to burs for fine interproximal reduction

For Removing Attachments:

Adhesive Removing Carbide Burs Strips residual composite off enamel after treatment

Polishing Burs Smooths and re-polishes enamel after composite removal

For Using Elastics:

Bondable Buttons | 10pk Bonded anchor points for attaching elastics to teeth

Bondable Buttons w/ Hooks | 10 pk Same as above, with an integrated hook for easier elastic placement

Elastics (recommended):

3/16, 2 ½ oz Light Elastics | 50 pk Smaller-diameter elastics for localized or tighter bite correction

1/4, 2 ½ oz Light Elastics | 50 pk Larger-diameter elastics for Class II/III correction with sustained light force

Other:

Pull tools Helps patients remove tight-fitting aligners without breaking them

Chewies Bite cylinders that help patients fully seat their aligners for accurate tracking

Cases Included in aligner delivery, extras on-hand are helpful for replacements
